

Lynx Mobile®

Enhancements – April 5th, 2013

On the evening of April 5th, 2013 McKesson Specialty Health will release Lynx Mobile Version 5.19. We are pleased to inform you that these enhancements, which are based on user feedback, offer the following new features and functionality:

Inventory Management Updates:

- **Item Attribute Settings for New Inventory** - Save time with our updated Item Attribute settings that flow to all locations when you initially add an item to inventory.
- **Internal Storage Location Item Transfers** – Easily track and transfer inventory between storage locations in your practice.

EMR Orders Integration Updates:

- **Improved Dispense & Queue visibility for multiple dispenses for the same medication** – View multiple dispenses for the same item separately when multiple order IDs or routes are indicated in the EMR Order.
- **Added Queue link** – When dispensing items within the EMR Orders tab, a [Queue] link appears next to the Item Name when it is already Queued for the same date of service, allowing you to navigate directly to the Queue.
- **Multiple Dispenses of the same Item in Dispense History** – Easily edit and view multiple dispenses of the same item.

Queue Updates:

- **Automated Queue Cleanup** – Configurable to suit your practice's needs, Queue entries will automatically be deleted when they are seven days beyond the date of service, providing cleaner data across reports that reference Queue data.

Reporting Updates:

- **Queue Report** - Now enhanced to include non-patient transactions, the Queue Report provides you with additional information about your on hand inventory levels, crucial to determining what you need to purchase in the coming days.
- **Refill Report** – Items can be sorted by storage location, allowing you to assess refill needs for multiple storage locations containing the same item.
- **Patient Data Quality Report** – Enhanced to display patients with invalid entries for height and weight.

Equipment Updates:

- **Remote Locking Device Alarm** - Now includes an option to customize and configure your remote locking device alarm settings. This will prevent the alarm from sounding during routine actions (Restock, Verify Inventory) while continuing to alert you to a critical situation when the lock is not properly latched.

Additional Improvements:

We've included a number of minor improvements and bug fixes that will improve your overall Lynx Mobile experience.

We're here to help

If you need assistance with Lynx Mobile, please call 800.482.6700 or email us at msh.techsupport@mckesson.com.

INVENTORY MANAGEMENT

Item Attribute Settings for New Inventory Items

While adding new items to inventory, item attributes (highlighted below in orange) will need to be selected only once, upon initial entry. The attributes selected will then be applied to the item at each storage location and/or across the multiple sites which you have access to, significantly reducing repetitive data entry as you add new items to your inventory.

Item Information

Item Name: * G-Csf 300 mcg Syringe 10 per Box
 Brand Name: * Neupogen 300 mcg Syringe 10 per Box

Item Code: * 55513-0924-10
 Inventory Unit: * Syringe
 Disp Units/Inv Unit: 300
 NDC/GCN: 55513-0924-10
 (based on item code) 45996

Catalog #: * 221-200
 Dispensing Unit: * mcg
 Billing Code: 11440

Report Seq No: 3
 Cost:
 Billing Unit: 300
 Billing Unit Override:

Active: Autowaste: Verify Count before Dispense: Allow Auto-Refill:
 Billable: Track Lot#: Witness on Dispense:

Inventory Information

Site	Storage Location	Auto-Refill	Label	Inv Qty	Partial Qty	Min	Max	Exp Date	DEL?
001	...Non Chemo Fridge	<input type="checkbox"/>	0	16	1	5	10		<input checked="" type="checkbox"/>
002	...Large Refrigerator	<input type="checkbox"/>	1	20	0	0	1		<input checked="" type="checkbox"/>
002	...Growth Factors	<input type="checkbox"/>	1	18	0	10	30		<input checked="" type="checkbox"/>
WEBUSOLM	...Refrigerator	<input type="checkbox"/>	1	51	250	0	1		<input checked="" type="checkbox"/>


Internal Storage Location Item Transfers

This release introduces internal storage location transfer capabilities so you can easily track and sort inventory amongst several storage locations within your practice. Storage locations are now available in the Transfer Location drop down list as destinations. (NOTE: You will need to add the item to the location first if the item has not been stored in the location you wish to transfer to) These internal transfers will be listed in the Queue for completion.

To complete a transfer:

- Navigate to Inventory Admin → select Inventory List
- Search for the item you wish to transfer
- Enter the amount to be transferred in the Transfer Amount field in the row corresponding to the Storage Location you wish to transfer from
- The Transfer Location drop down list will load

NOTE: If the location that you wish to transfer to does not appear, add the item to the destination location prior to conducting the transfer.

- Select the Storage Location you wish to transfer the Item to
- Click 

Inventory Administration

Generic Name: Item Code:

Brand Name: Catalog #:

Storage Location: Search Type:

[Inventory List](#) [Item Master](#) [Global Inventory Management](#) [Restock](#)

Prev 1 Next Note: You must save your changes before going to another page

Name Brand Name	Item Code Cat #	*Total Inv Total Min Total Max	Location	Inv Qty Adjustment Reason	Partial Qty Reason	Min Max	Restock Amount Restock Reason	Waste Amount Waste Reason	Auto Refill Exp. Date	Label	Transfer Amount	Transfer Location	DEL?
Dexamethasone 10 mg/1 mL SDV Dexamethasone 10 mg (25 vial/flat)	00641-0367- 25 S03-353 [New Location]	1988 Vial 0 Vial A 2 VialDrawer [change]	1988 Vial	0 mg	0 1	Vial	mg		1	25 Vial	<input type="text" value="Select One"/>	<input type="button" value="X"/>

* Refill amounts will be calculated from these values

Internal Locations
Cabinet 1
Door 1
Door 2
Drawer B
Drawer C
Drawer D
Door 3
Refrigerator
Procedure Codes
Fluids and Supplies
Short Inf Room
Drawer A
Drawer B
Drawer C
Drawer D
Door 1
Refrigerator
Procedure Codes
Fluids and Supplies

Copyright © 2013 McKesson Specialty Health. All rights reserved.

- Navigate to the Queue to and select the [Transfer](#) that you would like to execute
- Check the Dispense box → click **TRANSFER** to open the Storage Location to execute the transfer and complete the transfer by closing the Location or by clicking **CONFIRM** if required.

NOTE: The Queue will contain two entries for the internal Transfer – one to remove the Item from the Transfer From Location and one to restock the Item in the Transfer To Location

Queue - Select a set of transactions below to dispense, edit, delete, or view.

Patient			Date of Service		Transaction Type [all]			Other		
<input type="text"/>	<input type="text"/>	<input type="text"/>	Specific: 02/28/2013	-or- Range: <input type="text"/>	<input checked="" type="checkbox"/> Dispense	<input checked="" type="checkbox"/> Waste	<input type="checkbox"/> Only My Transactions			
<small>last</small>	<small>first</small>	<small>id</small>	<small>mm/dd/yyyy</small>		<input checked="" type="checkbox"/> Restock	<input checked="" type="checkbox"/> Transfer	<input checked="" type="radio"/> Patient View			
					<input type="radio"/> Location View					
					[reset all filters]					

Tx TYPE	Tx REASON	PATIENT NAME	USER ID	DATE OF SERVICE	LAST QUEUED TIME
Transfer	To Drawer D		aduprey		02/28/2013 12:01 PM
Transfer	From Drawer C		aduprey		02/28/2013 12:01 PM

Dispense:	Location	Item	Dispense Unit	Rt	Total In Queue	Inv Qty	Partial	Storage Location Override
<input checked="" type="checkbox"/>	Cabinet xyz	Drawer C Dexamethasone 4 mg per mL 30 mL Mdv - 00517-4930-25 (00517-4930-25)	2 outgoing	Vial	incoming 2,868	9 Vial	41 mg	

EMR ORDERS INTEGRATION IMPROVEMENTS

Dispensing Multiple Orders for the Same Medication

Lynx Mobile allows you to view multiple dispenses and queue transactions for the same item dispensed separately based on the order ID and route dispensed. This feature is now available on the Dispense-Add Items Page when EMR Orders Integration is enabled. The Dose-Dispense validation alert will appear when the dispense amount is not exactly equal to the ordered dose.

Dispense - Add Items [RETURN TO SEARCH](#)

Select items from the History, Treatment Regimen and/or Inventory menu.

Patient: Name not shown for patient confidentiality

DOB: 08/14/1958 **Height:** **Diagnosis:**
ID: 207129 **Weight:**
Physician: zzNeilan, Rosemary
Allergies:
Notes:

[History](#) [Treatment Regimen](#) [Items](#) **[EMR Orders](#)**

Order Date: 03/18/2013 [\[refresh\]](#) **ADD ALL** **Ordering Physician:** UNSPECIFIED

EMR Order Item Name	Ordered Dose	Order UOM	Lynx Mobile Item Name Brand or Code	Inv Qty	Partial
Hydrocortisone sod succinate, inj	100 mg		Hydrocortisone Sod Succ 100MG AOV SDV PF 1/EA Solu-Cortef 100MG AOV SDV PF 1/EA	5 Vial	0 mg
Fluorouracil, inj	480 mg		Fluorouracil 50 mg per mL 100 mL Fluorouracil 5000 mg SDV	996 Vial	3210 mg
Fluorouracil CIV, inj	2890 mg		Fluorouracil 50 mg per mL 100 mL Fluorouracil 5000 mg SDV	996 Vial	3210 mg
Diphenhydramine hcl, inj	25 mg		DiphenhydrAMINE 50 mg per mL 1 mL 00641-0376-25	22 Vial	30 mg

DEL ALL **Date of Service** **Clinical Trial** [QUEUE](#)

[PREV](#) 03/18/2013 [NEXT](#) NONE [NEW](#) [SAVE TO QUEUE](#)

Name (Brand or Code)	Route	Disp Amt	Inv Amt	Waste	Total	Inv Qty	Partial
Fluorouracil 50 mg per mL 100 mL - Fluorouracil 5000 mg SDV (10139-0063-01) [Instructions]	IV PUSH	480 mg		0 mg	480 mg	996 Vial	3210 mg
Fluorouracil 50 mg per mL 100 mL - Fluorouracil 5000 mg SDV (10139-0063-01) [Instructions]	CONTINUOUS INFUSION	2890 mg		0 mg	2890 mg	996 Vial	3210 mg

Viewing the Same Item Dispensed from Multiple Orders in the Queue

When EMR Orders Integration is enabled, dispense amounts of the same item will appear separately in the Queue as they relate to Order IDs or route, providing you with complete visibility into your dispensing activities.

Queue - Select a set of transactions below to dispense, edit, delete, or view.

Patient	Date of Service	Transaction Type [all]	Other																																													
EMR Orders Patient: Name not shown for patient confidentiality DOB: 01/01/1980 ID: 203595 Order Date: 11/05/2012 [refresh]																																																
<input checked="" type="checkbox"/> Waste <input type="checkbox"/> Only My Transactions <input checked="" type="checkbox"/> Transfer <input checked="" type="radio"/> Patient View <input type="radio"/> Location View																																																
DATE OF SERVICE ▼ LAST QUEUED TIME 11/05/2012 02/21/2013 01:54 PM																																																
<table border="1"> <thead> <tr> <th>EMR Order Item Name</th> <th>Physician</th> <th>Ordered Dose</th> <th>UOM</th> <th>Route</th> </tr> </thead> <tbody> <tr> <td> Lorazepam, inj</td> <td>Dunning,David M</td> <td>1 mg</td> <td></td> <td>IV PUSH</td> </tr> <tr> <td> Mesna, inj</td> <td>Dunning,David M</td> <td>1110 mg</td> <td></td> <td>I.V.</td> </tr> <tr> <td> Ifosfamide, inj</td> <td>Dunning,David M</td> <td>3324 mg</td> <td></td> <td>I.V.</td> </tr> <tr> <td> Mesna, inj</td> <td>Dunning,David M</td> <td>1110 mg</td> <td></td> <td>I.V.</td> </tr> <tr> <td> Mesna, inj</td> <td>Dunning,David M</td> <td>1110 mg</td> <td></td> <td>I.V.</td> </tr> <tr> <td> Dexamethasone sod phosphate, inj</td> <td>Dunning,David M</td> <td>10 mg</td> <td></td> <td>IV PUSH</td> </tr> <tr> <td> Etoposide, inj</td> <td>Dunning,David M</td> <td>180 mg</td> <td></td> <td>IV</td> </tr> <tr> <td> Lasix, inj (furosemide)</td> <td>Dunning,David M</td> <td>20 mg</td> <td></td> <td>IV PUSH</td> </tr> </tbody> </table>				EMR Order Item Name	Physician	Ordered Dose	UOM	Route	Lorazepam, inj	Dunning,David M	1 mg		IV PUSH	Mesna, inj	Dunning,David M	1110 mg		I.V.	Ifosfamide, inj	Dunning,David M	3324 mg		I.V.	Mesna, inj	Dunning,David M	1110 mg		I.V.	Mesna, inj	Dunning,David M	1110 mg		I.V.	Dexamethasone sod phosphate, inj	Dunning,David M	10 mg		IV PUSH	Etoposide, inj	Dunning,David M	180 mg		IV	Lasix, inj (furosemide)	Dunning,David M	20 mg		IV PUSH
EMR Order Item Name	Physician	Ordered Dose	UOM	Route																																												
Lorazepam, inj	Dunning,David M	1 mg		IV PUSH																																												
Mesna, inj	Dunning,David M	1110 mg		I.V.																																												
Ifosfamide, inj	Dunning,David M	3324 mg		I.V.																																												
Mesna, inj	Dunning,David M	1110 mg		I.V.																																												
Mesna, inj	Dunning,David M	1110 mg		I.V.																																												
Dexamethasone sod phosphate, inj	Dunning,David M	10 mg		IV PUSH																																												
Etoposide, inj	Dunning,David M	180 mg		IV																																												
Lasix, inj (furosemide)	Dunning,David M	20 mg		IV PUSH																																												
<table border="1"> <thead> <tr> <th>Dispense? Location</th> <th>Total in Queue</th> <th>Inv Qty</th> <th>Partial</th> <th>Storage Location Override</th> <th>Instructions</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/> Cabinet 1</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td><input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)</td> <td>1,110 dispensed</td> <td>mg</td> <td>I.V.</td> <td>3,330</td> <td>20 Vial 0 mg [Instructions]</td> </tr> <tr> <td><input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)</td> <td>1,110 dispensed</td> <td>mg</td> <td>I.V.</td> <td>3,330</td> <td>20 Vial 0 mg [Instructions]</td> </tr> <tr> <td><input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)</td> <td>1,110 dispensed</td> <td>mg</td> <td>I.V.</td> <td>3,330</td> <td>20 Vial 0 mg [Instructions]</td> </tr> </tbody> </table>				Dispense? Location	Total in Queue	Inv Qty	Partial	Storage Location Override	Instructions	<input checked="" type="checkbox"/> Cabinet 1						<input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)	1,110 dispensed	mg	I.V.	3,330	20 Vial 0 mg [Instructions]	<input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)	1,110 dispensed	mg	I.V.	3,330	20 Vial 0 mg [Instructions]	<input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)	1,110 dispensed	mg	I.V.	3,330	20 Vial 0 mg [Instructions]															
Dispense? Location	Total in Queue	Inv Qty	Partial	Storage Location Override	Instructions																																											
<input checked="" type="checkbox"/> Cabinet 1																																																
<input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)	1,110 dispensed	mg	I.V.	3,330	20 Vial 0 mg [Instructions]																																											
<input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)	1,110 dispensed	mg	I.V.	3,330	20 Vial 0 mg [Instructions]																																											
<input checked="" type="checkbox"/> Door 1 Mesna 1000 mg/10 mL MDV - Mesna 100 mg per mL 10 mL Mdv (63323-0733-10)	1,110 dispensed	mg	I.V.	3,330	20 Vial 0 mg [Instructions]																																											

Viewing Multiple Dispenses of the Same Item from the Dispense History

Lynx Mobile now allows you to view multiple dispenses of the same item separately in the Dispense History so you can easily view and edit your dispenses. As an added feature to simplify your experience, the icon next to the Name (Brand or Code) will allow you to view the Order ID of the original EMR order used for the Dispense.

A [\[Queue\]](#) link appears when an Item is Queued for the same date of service.

When dispensing items within the EMR Orders tab, a [\[Queue\]](#) link appears next to the Item Name when it is already Queued for the same date of service, allowing you to navigate directly to the Queue.

Clicking on the [\[Queue\]](#) link will allow you to quickly navigate directly to the Queue page to view the patient's Dispense details for the Date of Service in focus.

Dispense - Add Items RETURN TO SEARCH

Select items from the History, Treatment Regimen and/or Inventory menu.

Patient: Name not shown for patient confidentiality

DOB: 07/12/1961 Height: 70 in Diagnosis:
 ID: 236986 Weight: 199 lbs
 Physician: zzNeilan, Rosemary
Allergies:
Notes:

[History](#) [Treatment Regimen](#) [Items](#) [EMR Orders](#)

Order Date: 12/05/2012 [\[refresh\]](#) + ADD ALL Ordering Physician: UNSPECIFIED

EMR Order Item Name	Ordered Dose	Order UOM	Lynx Mobile Item Name Brand or Code	Inv Qty	Partial
Pamidronate disodium, inj	90	mg	Pamidronate Disodium 9 mg per mL 10 mL 90 mg SDV 61703-0326-18 [Queue]	6 Vial	0 mg
Hydrocortisone sod succinate, inj	100	mg	Hydrocortisone Sod Succ 100MG AOV SDV PF 1/EA Solu-Cortef 100MG AOV SDV PF 1/EA [Queue]	5 Vial	0 mg
Methylprednisolone sod succ, inj	125	mg	MethylPREDNISolone SOD SUCC Solu-Medrol 1 gm AOV SDV [Queue]	7 Vial	0.5 G
Dinhenhydramine hdl. inj	25	mg	DinhenhydrAMINE 50 mg per mL 1 mL	22 Vial	30 mg

DEL ALL **Date of Service** **Clinical Trial** QUEUE

← PREV
12/05/2012
NEXT →
NONE
NEW
SAVE TO QUEUE

Name (Brand or Code)	Route	Disp Amt	Inv Amt	Waste	Total	Inv Qty	Partial

QUEUE

Queued Entries older than 7 days past the date of service will be cleaned up nightly

Configurable to suit your practice's needs, Queue entries will automatically be deleted when they are seven days beyond the date of service, providing cleaner data across reports that reference Queue data.

To adjust the days that your queue entries will be saved, go to Admin → Practice Info and scroll down to the Queue Settings section to adjust the Save Queue Entries.

Practice Info

This page is displaying the details specific to your practice.

General Info

Practice Name: * Lynx Mobile Oncology

Addr 1: * 4565 McKesson St

Addr 2:

City: * San Francisco

State: * CA Zip: * 94104

Time Zone: US/Central

Phone: 1-800-482-6700

FAX: 800-482-6700

Contact Name: Administrator

Contact e-mail: Administrator@practice.com

Tax ID: 12345678

Auto Logoff: * 100 minutes Passwords Expire After: * 365 days

Force Inventory Count Before Dispense:

Printed Label Settings

Display Transaction ID: Display Insurance:

Queue Settings

Queue View: Patient

Default Date of Service: Today+Past

Save Queue Entries: * 7 days past their DOS
Min 0, Max 30

REPORTING

Queue Report Enhanced to Include Non Patient Transactions

Now enhanced to include non-patient transactions, the Queue Report provides you with additional information about your on hand inventory levels and demand, crucial to determining what you need to purchase in the coming days. The non-patient transactions will include any Transfers, Restocks and Waste transactions.

NOTE: Please be advised that a refill item sent to the Queue to restock should be restocked from the Queue prior to running the Queue Report; otherwise the restock amount will not be deducted from the Placed Order Total.

Lynx Mobile - Queue Report (Without Patients)

Downloaded 03/28/2013 11:33 AM (PDT)

Queue with Date of Service Between **03/01/2013** and **04/01/2013**

DiphenhydrAMINE 50 mg per mL 1 mL 00641-0376-25 (00641-0376-25)		
Site	Service Date	Disp. Amt Vial
ARCATACALM	03/28/13	10
	Queued Patient Total	10
	On Hand Total	2.1
	+ Queued Outbound Transfers	5
	- Queued Restocks	5
	+ Queued Inventory Waste	0.1
	Additional Amt Needed	8

Refill Report

Items can be sorted by storage location, allowing you to assess refill needs for multiple storage locations containing the same item.

When the Consolidate Locations checkbox is left unchecked, the option to sort by Storage Location will be available for items in the Refill Report. Simply click [Storage Location](#) to sort by location.

Refill Report
The Refill Report displays the list of items in your inventory that are at or below the minimum inventory level.

Report Parameters

Show Auto-Refill Items
 Show Non-Auto-Refill Items
 Consolidate Locations:

Report Results

[PRINT DETAILS](#) [DOWNLOAD EXCEL](#) [RESET STATUS](#)

Name Brand Name	Item Code Cat #	Storage Location	Min Max	Inv Qty	Status Order #	Requested Date Requested Qty	Ordered Date Ordered Qty	Restocked Date Restocked Inv Qty	Warehouse Availability	Reset?
Bortezomib 3.5MG Powder Velcade 3.5MG Powder	63020-0049-01 910-035	Cabinet xyz - Door 1	70 Vial 100 Vial	5 Vial	Refill Requested	03/04/2013 02:10 PM CST 33 Vial			D	<input type="checkbox"/>
Darbepoetin Alfa 300 mcg Vial Aranesp 300 mcg per 1 mL SDV Albumin- Free	55513-0110-01 226-301	Chemo Fridge	0 Vial 1 Vial	0 Vial	Not Yet Ordered				✓	<input type="checkbox"/>
Darbepoetin Alfa 300 mcg Vial Aranesp 300 mcg per 1 mL SDV Albumin- Free	55513-0110-01 226-301	daRefrigerator	0 Vial 1 Vial	0 Vial	Not Yet Ordered				✓	<input type="checkbox"/>
LORazepam 2 mg per mL Sol 1 mL Fliptop vial LORazepam 2 mg per mL Sol 1 mL Fliptop vial	00409-6778-02 261-104	Non Chemo Fridge1	4 Vial 8 Vial	0 Vial	Refill Requested	03/04/2013 02:10 PM CST 8 Vial			✓	<input type="checkbox"/>
MitoMYCIN 20mg SDV PWD MitoMYCIN 20mg SDV PWD	16729-0108-11 649-891	Rheumatology Clinic	0 vial 1 vial	0 vial	Refill Requested	02/26/2013 02:10 PM CST 1 vial			✓	<input type="checkbox"/>
Acetaminophen 325 mg Tablet 100/Bottle Acetaminophen 325 mg Tablet 100/Bottle	00536-3222-01 350-421		10 Bottle 20 Bottle	8 Bottle	Refill Requested	02/26/2013 02:10 PM CST 12 Bottle			Unspecified	<input type="checkbox"/>

[PRINT DETAILS](#) [DOWNLOAD EXCEL](#) [RESET STATUS](#)

Patient Data Quality Report Data Validation

Patients with invalid entries for Height and Weight will appear on the worklist until they are updated with appropriate values and saved.

Customize and configure your remote locking device alarm settings

You can prevent your remote locking device alarm from sounding during routine activities, such as restocking or verifying inventory, by customizing your alarm settings. The current default is 2 minutes; however we added a feature that allows you to configure the amount of time to suit your practice's preference.

- To access the Storage Location Configuration click on Admn → Storage Location Configuration and select the location that you would like to customize.
- Enter the delay time in seconds and click **SAVE**

Storage Location Configuration

To add a storage location, select a storage type and a storage name from the tree view, and click Add button.
 To delete a storage location, select a storage name from the tree view and click Delete button.
 To move a storage location, select a storage name from the tree view and click arrow up or down button.
 To save the changes click Save button.
 Cabinets Attached: (none)

Storage Types

- Full Cabinet
- Half Cabinet
- Cabinet
- Door
- Shelf
- Drawer
- Pocket
- RemoteLock
- Unsecured

ADD DELETE

↑ iKnowMed East Bay Oncology 6.4

- Cabinet 1
 - Door 1
 - Door 2
 - Drawer A
 - Drawer B
 - Drawer C
 - Drawer D
 - Door 3
- Refrigerator
- Procedure Codes
- Fluids and Supplies

↑

↓

Name:

Cabinet #:

Alarm Delay: seconds

⚙ Min 0, Max 240

⚙ Leaving this value blank causes the hardware default to be used

⚙ Not all refrigerator locks respect this setting

Require Hardware for Dispenses:

SAVE CANCEL

Additional Improvements

- **Label Instructions:** Added or modified Label Instructions will be saved if you navigate away from the Queue.
- **Queue Inventory Qty:** The Queue will now show the actual amount on hand, in units dispensed, when there is insufficient Inventory to complete the transaction. The out of stock message will indicate the insufficient amount remaining across multiple Storage Locations (if applicable).

Queue - Select a set of transactions below to dispense, edit, delete, or view.

Tx TYPE	Tx REASON	PATIENT NAME	USER ID	DATE OF SERVICE	LAST QUEUED TIME
Dispense		Anderson, Nancy June (DOB: 09/18/1937, ID: 211059)	LAGLASS	02/20/2013	02/20/2013 05:26 PM

Patient: Anderson, Nancy June [More Patient Info / BSA Calculator](#) [View EMR Orders](#) Date of Service: 02/20/2013
 DOB: 09/18/1937 Height: 58 in Supervising Physician: Bufill, Jose A (6)
 ID: 211059 Weight: 95 lbs Clinical Trial: NONE
 Physician: Bufill, Jose A Diagnosis: 204.00 - ACT LYM LEUK W/O RMSION, Stage: FAB-L2
 Allergies: Amoxicillin, Penicillins, erythromycin
 Notes:

Dispense?	Location	Item	Dispense Unit (Waste/Total)	Rt	Total In Queue	Inv Qty	Partial	Lot # Exp Date	Storage Location Override Instructions
<input type="checkbox"/>	Cabinet xyz	Drawer C CARBOPLATIN 10MG/ML 60ML MDV -	600 dispensed mg		600	0 Vial	450 mg		[Instructions]

Only 450 mg available across 2 storage locations

- **Browsers:** Lynx Mobile is now optimized for use with Chrome and Firefox browsers.
- **Queue Search:** Search performance in the Queue, when searching by Location View, is now optimized.
- **Patient Search:** We improved patient search performance when more than 100 name matches are found.
- **SuperBill Report:** The SuperBill will now reflect the total amount of an item dispensed (regardless of the route used) rather than itemizing each Dispense by route, to ease the calculation of billable units for the day.
- **Queue Report:** Queue Report items are now organized by Item Code rather than Item Name to prevent combining distinct items with identical names.
- **Global Label Item Association:** The Global Label Item Association page will now include a scrolling pane to handle long lists.